


# MEDICAL PACKAGING

## BOLSAPACK • PACKAGING FOR STEAM/EO/FORMOL

Pouches and rolls, manufactured with 60 g/m<sup>2</sup> medical paper and a plastic film (blue coloured). Indicators for Steam/EO/Formol are also printed in packaging.

### PAPER/FILM FLAT BAG

Code	Sizes		U/Box	Box/cm	Box/Kg
MBO1025	75 x 150 mm	3,0" x 5,9"	2.700	15 x 27 x 15	3,5
MBO1000	75 x 210 mm	3,0" x 8,3"	2.500	31 x 29 x 11	4,8
MBO1001	75 x 280 mm	3,0" x 11,0"	2.400	31 x 29 x 11	6,5
MBO1003	100 x 200 mm	3,9" x 7,9"	2.400	41 x 29 x 11	5,9
MBO1004	100 x 240 mm	3,9" x 9,4"	2.400	41 x 29 x 11	7
MBO1005	100 x 280 mm	3,9" x 11,0"	2.400	41 x 29 x 11	8,5
MBO1026	100 x 300 mm	3,9" x 11,8"	2.500	41 x 43 x 11	9
MBO1006	100 x 420 mm	3,9" x 16,5"	2.400	41 x 43 x 11	12,5
MBO1007	150 x 210 mm	5,9" x 8,3"	2.000	32 x 43 x 11	8
MBO1008	150 x 280 mm	5,9" x 11,0"	2.000	31 x 57 x 11	10,5
MBO1009	150 x 320 mm	5,9" x 12,6"	2.000	37 x 65 x 11	11,5
MBO1010	150 x 360 mm	5,9" x 14,2"	2.000	37 x 65 x 11	13
MBO1021	150 x 400 mm	5,9" x 15,7"	1.200	32 x 43 x 11	8,7
MBO1011	200 x 320 mm	7,9" x 12,6"	1.800	37 x 65 x 11	14
MBO1012	200 x 360 mm	7,9" x 14,2"	1.800	37 x 65 x 11	15,8
MBO1013	200 x 420 mm	7,9" x 16,5"	1.200	41 x 43 x 11	12,5
MBO1023	200 x 500 mm	7,9" x 19,7"	1.200	51 x 51 x 11	14,8
MBO1014	250 x 360 mm	9,8" x 14,2"	1.200	51 x 43 x 11	13,5
MBO1015	250 x 420 mm	9,8" x 16,5"	1.200	51 x 43 x 11	15,5
MBO1016	250 x 500 mm	9,8" x 19,7"	1.200	51 x 51 x 11	18,5
MBO1017	300 x 420 mm	11,8" x 16,5"	600	32 x 43 x 11	9,5
MBO1018	300 x 500 mm	11,8" x 19,7"	600	31 x 57 x 11	11
MBO1024	300 x 580 mm	11,8" x 22,8"	700	41 x 59 x 11	14,2
MBO1019	400 x 580 mm	15,7" x 22,8"	600	41 x 59 x 11	17



### PAPER/FILM GUSSET BAG

Code	Sizes		U/Box	Box/cm	Box/Kg
MBF1202	200 x 55 x 400 mm	7,9 x 2,2" x 15,7"	600	43 x 26 x 23	7

# MEDICAL PACKAGING

## BOLSAPACK • PACKAGING FOR STEAM/EO/FORMOL

Pouches and rolls, manufactured with 60 g/m<sup>2</sup> medical paper and a plastic film (blue coloured). Indicators for Steam/EO/Formol are also printed in packaging.

### PAPER/FILM SELFSEALING BAG

Code	Sizes		U/Pack	U/Box	Box/cm	Box/Kg
MBA1080	75 x 200 mm	3,0" x 7,9"	200	1.200	15 x 24 x 17	2,7
MBA1081	90 x 200 mm	3,5" x 7,9"	200	1.200	28 x 27 x 12	4
MBA1082	90 x 230 mm	3,5" x 9,1"	200	1.200	28 x 27 x 12	4,6
MBA1084	100 x 240 mm	3,9" x 9,4"	200	1.200	31 x 29 x 11	4,5
MBA1085	130 x 250 mm	5,1" x 9,8"	200	1.200	41 x 29 x 11	5,8
MBA1086	130 x 330 mm	5,1" x 13,0"	200	1.200	41 x 37 x 11	7,6
MBA1087	140 x 280 mm	5,5" x 11,0"	200	1.200	43 x 32 x 11	7
MBA1088	200 x 330 mm	7,9" x 13,0"	200	800	41 x 37 x 11	7,8
MBA1089	250 x 400 mm	9,8" x 15,7"	200	800	51 x 43 x 11	12
MBA1090	300 x 390 mm	11,8" x 15,4"	***	400	32 x 43 x 11	6,9
MBA1091	300 x 450 mm	11,8" x 17,7"	***	400	31 x 49 x 11	8



### PAPER/FILM FLAT ROLL

Code	Sizes		U/Box	Box/cm	Box/Kg
MRP1040	50 mm x 200 m	2,0" x 656 Feet	6	31 x 25 x 22,5	6
MRP1041	75 mm x 200 m	3,0" x 656 Feet	4	31 x 25 x 22,5	8
MRP1042	100 mm x 200 m	3,9" x 656 Feet	4	43 x 26 x 22,5	11
MRP1043	125 mm x 200 m	4,9" x 656 Feet	4	43,5 x 32 x 22,5	14
MRP1044	150 mm x 200 m	5,9" x 656 Feet	4	43,5 x 32 x 22,5	17
MRP1045	200 mm x 200 m	7,9" x 656 Feet	2	43 x 26 x 22,5	11
MRP1046	250 mm x 200 m	9,8" x 656 Feet	2	52 x 2 8x 22,5	14
MRP1047	300 mm x 200 m	11,8" x 656 Feet	2	43,5 x 32 x 22,5	17
MRP1050	350 mm x 200 m	13,8" x 656 Feet	1	43 x 26 x 22,5	10
MRP1048	400 mm x 200 m	15,7" x 656 Feet	1	43 x 26 x 22,5	11
MRP10491	500 mm x 100 m	19,7" x 328 Feet	1	53 x 18 x 18	7

### PAPER/FILM GUSSET ROLL

Code	Sizes		U/Box	Box/cm	Box/Kg
MRF1059	75 x 25 mm x 100 m	3,9 x 1,6" x 328 Feet	4	31 x 25 x 22,5	6
MRF1060	100 x 40 mm x 100 m	3,9 x 1,6" x 328 Feet	4	43 x 26 x 22,5	7
MRF1061	150 x 40 mm x 100 m	5,9 x 1,6" x 328 Feet	4	43,5 x 32 x 22,5	11
MRF1062	200 x 60 mm x 100 m	7,8 x 2,3" x 328 Feet	2	43 x 26 x 22,5	7
MRF1063	250 x 60 mm x 100 m	9,8 x 2,3" x 328 Feet	2	52 x 28 x 22,5	8
MRF1064	300 x 80 mm x 100 m	11,8 x 3,1" x 328 Feet	1	31 x 25 x 22,5	6
MRF1065	400 x 80 mm x 100 m	15,7 x 3,1" x 328 Feet	1	43 x 26 x 22,5	7