

MEDICAL PACKAGING

BOLSAPACK • PACKAGING FOR LOW TEMPERATURE (VH₂O₂) - EO

Pouches and Rolls manufactured with Tyvek® and transparent film. Printed with chemical indicator for low temperature (Peroxide - VH₂O₂).

TYVEK®/FILM FLAT BAG

Code	Sizes		U/Box	Box/Cm	Box/Kg
MBT1160	75 x 210 mm	3 x 8"	1.000	16 x 22 x 15	2,4
MBT1161	100 x 280 mm	4 x 11"	1.000	16 x 22 x 15	4,1
MBT1162	150 x 320 mm	6 x 13"	1.000	31 x 33 x 15	7,0
MBT1163	200 x 420 mm	8 x 17"	350	43 x 21 x 11	4,0
MBT1164	250 x 500 mm	10 x 20"	350	28 x 60 x 12	6,5
MBT1165	320 x 560 mm	13 x 22"	400	41 x 59 x 11	10,5
MBT1166	400 x 580 mm	16 x 23"	350	41 x 59 x 11	12,0

TYVEK®/FILM FLAT ROLL

Code	Sizes		U/Box	Box/Cm	Box/Kg
MRT1171	75 mm x 70 m	3" x 230 Feet	4	31 x 25 x 23	3,5
MRT1172	100 mm x 70 m	4" x 230 Feet	4	41 x 21 x 22	4,5
MRT1174	150 mm x 70 m	6" x 230 Feet	4	37 x 34 x 22	6,5
MRT1175	200 mm x 70 m	8" x 230 Feet	2	41 x 21 x 22	4,5
MRT1176	250 mm x 70 m	10" x 230 Feet	2	43 x 26 x 23	6,0
MRT1177	300 mm x 70 m	12" x 230 Feet	2	37 x 21 x 34	6,0
MRT1180	350 mm x 70 m	14" x 230 Feet	1	41 x 22 x 23	3,6
MRT1178	400 mm x 70 m	16" x 230 Feet	1	41 x 21 x 22	4,5
MRT1179	500 mm x 70 m	20" x 230 Feet	1	53 x 18 x 18	6,0




Tyvek®/Film bag



Opening instructions



Tyvek®/Film roll



Plasma indicator Tyvek®

TYVEK®/FILM SELFSEALING POUCH

Selfsealing bags manufactured with Tyvek® and transparent film. Printed with chemical indicator for low temperature (Peroxide - VH₂O₂) and chemical indicator for EO. Sealed with adhesive tape, no sealing machine required.

Code	Sizes		U/Box	U/Case	Box/Cm	Box/Kg
MBTA1180	75 x 200 mm	3 x 8"	900	150	17 x 25 x 17	2,5
MBTA1182	90 x 230 mm	4,5 x 9"	900	150	22 x 28 x 17	3,2
MBTA1184	100 x 240 mm	4 x 9,5"	900	150	23 x 29 x 17	3,7
MBTA1185	130 x 250 mm	5 x 10"	900	150	28 x 30 x 17	4,8
MBTA1187	140 x 280 mm	5,5 x 11"	900	150	30 x 33 x 17	5,7
MBTA1188	200 x 330 mm	8 x 13"	600	100	42 x 38 x 12	6,2
MBTA1189	250 x 400 mm	10 x 16"	600	100	52 x 45 x 12	9,2
MBTA1191	300 x 450 mm	12 x 18"	300	***	31 x 50 x 11	6,0



Tyvek®/Film selfsealing bag



With selfsealing tape



Instructions for use



2 indicators: VH₂O₂ and EO

Tyvek® is a registered trade mark of DUPONT