



High resistance bag with 3 sealings to ensure a good sealing strength. With thumbholes in both sides for an easy opening



Manufactured with soft and malleable materials, to easy fold it without any marks

REINFORCED POUCH

High resistance pouches, in blue colour for the packing of big and heavy products. This kind of bags offer a high tearing resistance and high permeability.

Manufactured with one side non-woven and one side film and conforms to ISO 11607-1 and EN 868-1.

It is considered a single use product and can be sterilized by Steam, Ethylene Oxide or Formaldehyde. Indicators for each process are printed in pouches.

Available in the following sizes:

Code	Sizes	U/Box	Box/Cm	Box/kg
MBCOVER001	200 x 420 mm 8" x 17"	1000	41 x 43 x 11	10,0
MBCOVER002	270 x 360 mm 11" x 14"	1000	28 x 38 x 17	11,5
MBCOVER003	270 x 450 mm 11" x 18"	500	47 x 57 x 14	7,5
MBCOVER004	300 x 500 mm 12" x 20"	500	31 x 57 x 11	9,5
MBCOVER005	400 x 520 mm 16" x 20"	500	41 x 59 x 11	12,5
MBCOVER006	400 x 580 mm 16" x 23"	500	41 x 59 x 11	14,5
MBCOVER007	500 x 600 mm 20" x 24"	500	51 x 61 x 13	18,5
MBCOVER008	570 x 720 mm 22" x 28"	250	41 x 59 x 23	14,0



Printed with 3 indicators for sterilization by STEAM, ETO and FORMALDEHYDE

CHARACTERISTICS

- High resistance
- High permeability
- Special for heavy and bulky products
- Wide range of sizes

bolsaplastmedical