



High resistance bag with 3 sealings to ensure a good sealing strength. With thumbholes in both sides for an easy opening



Manufactured with soft and malleable materials, to easy fold it without any marks

## REINFORCED POUCH

High resistance pouches, in blue colour for the packing of big and heavy products. This kind of bags offer a high tearing resistance and high permeability.

Manufactured with one side non-woven and one side film and conforms to ISO 11607-1 and EN 868-1.

It is considered a single use product and can be sterilized by Steam, Ethylene Oxide or Formaldehyde. Indicators for each process are printed in pouches.

Available in the following sizes:

| Code       | Sizes                  | U/Box | Box/Cm       | Box/kg |
|------------|------------------------|-------|--------------|--------|
| MBCOVER001 | 200 x 420 mm 8" x 17"  | 1200  | 41 x 59 x 23 | 13,0   |
| MBCOVER002 | 270 x 360 mm 11" x 14" | 500   | 31 x 51 x 23 | 6,4    |
| MBCOVER003 | 270 x 450 mm 11" x 18" | 500   | 47 x 57 x 14 | 7,5    |
| MBCOVER004 | 300 x 500 mm 12" x 20" | 500   | 31 x 57 x 11 | 9,5    |
| MBCOVER005 | 400 x 520 mm 16" x 20" | 500   | 41 x 59 x 11 | 12,5   |
| MBCOVER006 | 400 x 580 mm 16" x 23" | 500   | 41 x 59 x 11 | 14,5   |
| MBCOVER007 | 500 x 600 mm 20" x 24" | 500   | 51 x 61 x 13 | 18,5   |
| MBCOVER008 | 570 x 720 mm 22" x 28" | 250   | 41 x 59 x 23 | 14,0   |



Printed with 3 indicators for sterilization by STEAM, ETO and FORMALDEHYDE

## CHARACTERISTICS

- High resistance
- High permeability
- Special for heavy and bulky products
- Wide range of sizes

**bolsaplast**medical