

# BOLSADUST



Product information and instructions for use in several languages



Inviolable sealing, to make impossible reopening



Totally transparent, to see product inside without opening bag. Perfect view of inside packaging. A white square is printed, so user can write in

## CHARACTERISTICS

- Added protection to sterile pouches
- Excellent transparency
- Texts in several languages
- Self seal clousure
- With a white mark to write in
- Large flap for easy sealing
- Latex free
- Non-slip material
- Wide range of sizes

## PROTECTIVE BAG

Clear plastic bag designed to protect and to prevent deterioration of sterile packaging. It can be also used as an extra protection during transport of such kind of packaging, to avoid breaks and damages. This bag has printed in a white section where it is possible to write additional information. It is manufactured with non-slip PE and it is closed with an inviolable adhesive tape. It is considered a single use product and manufactured according ISO11607 standard.

This product is available in the following sizes:

Code	Sizes		U/Box	Box/cm	Box/Kg
MBDUST001	150 x 250 mm	6 x 10"	500	56 x 45 x 4	2,7
MBDUST002	200 x 300 mm	8 x 12"	500	56 x 45 x 4	4,3
MBDUST003	250 x 400 mm	10 x 16"	250	56 x 45 x 4	3,5
MBDUST004	300 x 500 mm	12 x 20"	250	61 x 45 x 4	5,4
MBDUST005	400 x 550 mm	16 x 22"	100	43 x 45 x 4	2,8
MBDUST006	400 x 760 mm	16 x 30"	100	43 x 45 x 4	3,8
MBDUST007	420 x 700 mm	17 x 28"	100	43 x 45 x 4	3,8
MBDUST008	500 x 670 mm	20 x 26"	100	56 x 45 x 4	4,8
MBDUST009	550 x 750 mm	22 x 30"	100	56 x 45 x 4	5,5
MBDUST010	600 x 760 mm	24 x 30"	100	61 x 45 x 4	5,9