


Product information and instructions for use in several languages



Inviolable sealing, to make impossible reopening

CHARACTERISTICS

- Added protection to sterile pouches
- Excellent transparency
- Texts in several languages
- Self seal clousure
- With a white mark to write in
- Large flap for easy sealing
- Latex free
- Non-slip material
- Wide range of sizes



Totally transparent, to see product inside without opening bag. Perfect view of inside packaging. A white square is printed, so user can write in

PROTECTIVE BAG

Clear plastic bag designed to protect and to prevent deterioration of sterile packaging. It can be also used as an extra protection during transport of such kind of packaging, to avoid breaks and damages. This bag has printed in a white section where it is possible to write additional information. It is manufactured with non-slip PE and it is closed with an inviolable adhesive tape. It is considered a single use product and manufactured according ISO11607 standard.

This product is available in the following sizes:

Code	Sizes		U/Box
MBDUST001	150 x 250 mm	5.9" x 9.8"	500
MBDUST002	200 x 300 mm	7.9" x 11.8"	500
MBDUST003	250 x 400 mm	9.8" x 15.7"	250
MBDUST004	300 x 500 mm	11.8" x 19.7"	250
MBDUST005	400 x 550 mm	15.7" x 21.6"	100
MBDUST006	400 x 760 mm	15.7" x 30.0"	100
MBDUST007	420 x 700 mm	16.5" x 27.5"	100
MBDUST008	500 x 670 mm	19.7" x 26.3"	100
MBDUST009	550 x 750 mm	21.6" x 29.5"	100
MBDUST010	600 x 760 mm	23.6" x 30.0"	100