

BOLSAPOCKET



Pockets perfectly located for an easy introduction of tools. Several window sizes to make it easy to use. Printed with indicators for Steam, ETO and Formaldehyde



Soft and malleable material, to easy fold it without any marks

SORTING BAG

Bags with pockets, in blue colour and transparent windows. Specially designed for sorting products (tools, instruments, etc.) inside the containers, flexible packaging (bags and rolls) or in the sterilization areas. They can also be used as an extra protection during instrumental transport, to avoid any damage or accidental breaking. By using the different pockets, materials can be easily classified and storage and distribution is optimized. It is considered a single use product and it is important to remark that it is not a sterile barrier system.


This product is available in the following size:

Code	Size	U/Box	Box/Kg
MBPOCKET001	500 x 210 mm 19.7" x 8.3"	300	3,7

CHARACTERISTICS

- Several sizes of pockets
- Easy visual identification of instrumental
- It can be rolled up to safe space

bolsaplastmedical



Film resistant to high temperatures, with wide sealings to avoid breaking when filling it. Windows are big enough to make easy introduction and removal of materials to be sterilized