

# BOLSACREPE



Excellent quality crepe paper in green colour. Soft and malleable material, to easy fold it without any marks.



For single use and for double wrapping.

## GREEN CREPED PAPER

Green creped paper sheets of 60 g for wrapping and sterilization of medical products. It is specially used in hospitals, medical centres, veterinarian centres, CSSD, etc.

It is also used as tray protector and for double wrapping. Manufactured from 100% cellulose and conforms to ISO 11607-1 and EN 868-2. Basically for wrapping medical devices to be sterilized.

It is considered a single use product and can be sterilized by Steam, Ethylene Oxide, Formaldehyde or Irradiation.

This product is available in the following sizes:

## CHARACTERISTICS

- Soft and malleable
- Excellent bacterial barrier
- Liquid repellent
- Wide range of sizes

**bolsaplast**medical

Code	Sizes		U/Box	Box/cm	Box/Kg
MCREPV500	500 x 500 mm	19 x 19"	750	52 x 52 x 11	11,2
MCREPV600	600 x 600 mm	23 x 23"	750	62 x 62 x 11	16,2
MCREPV750	750 x 750 mm	29 x 29"	348	82 x 41 x 11	11,8
MCREPV900	900 x 900 mm	35 x 35"	252	92 x 38 x 12	12,2
MCREPV1000	1000 x 1000 mm	39 x 39"	252	102 x 38 x 12	15,0
MCREPV1200	1200 x 1200 mm	47 x 47"	126	122 x 38 x 85	10,0